



**Product Specification**  
**Multi-purpose container, 26 ml, (LxØ): 45 x 70 mm, PP, transparent**

**Page 1**



**Product description**

Order number	75.568
Product description	Multi-purpose container, max. working volume: 26 ml, (LxØ): 45 x 70 mm, Ø opening: 70 mm, transparent, material: PP

**Product characteristics**

Graduation	no
Application	collection and storage of urine
Volume of work	26 ml

**Size**

Sample volume	26 ml
Max. volume	26 ml
Diameter	70 mm
Diameter of opening	70 mm
Height of product	45 mm
Length of product	45 mm
Length excluding cap	45 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Multi-purpose container, 26 ml, (LxØ): 45 x 70 mm, PP, transparent**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0036 kg
Weight of case	4.25 kg
EAN of inner box	4038917082549
EAN case	4038917001625

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

