




Product Specification
**Snap-on lid, without spout, PE, red, 1,000
piece(s)/bag**

Page 1



Product description

Order number	76.1414
Product description	Snap-on lid, without spout, suitable for 125-ml container 75.1355, material: PE, red, 1,000 piece(s)/bag

Product characteristics

Closure type	snap-on lid
Design	without spout
Suitable for	125-ml container 75.1355

Material & colours

Product material	Polyethylene (PE)
Colour of product	red

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Snap-on lid, without spout, PE, red, 1,000
piece(s)/bag

Page 2

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of box	380 mm
Width of box	190 mm
Height of box	195 mm
Depth of case	590 mm
Width of case	390 mm
Height of case	285 mm
Case volume	0.0656 cbm
Weight of product	0.0024 kg
Weight of case	3.14 kg
EAN of inner box	4038917086554
EAN case	4042836001499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT