



**Product Specification**  
**Snap-on lid, with spout, PE, red, 1,000**  
**piece(s)/bag**

**Page 1**



**Product description**

Order number	76.1477
Product description	Snap-on lid, with spout, suitable for 125-ml container 75.1355, material: PE, red, 1,000 piece(s)/bag

**Product characteristics**

Closure type	snap-on lid
Design	with spout
Suitable for	125-ml container 75.1355

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	red

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Snap-on lid, with spout, PE, red, 1,000**  
**piece(s)/bag**

**Page 2**

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of box	420 mm
Width of box	320 mm
Height of box	140 mm
Depth of case	590 mm
Width of case	390 mm
Height of case	265 mm
Case volume	0.061 cbm
Weight of product	0.0037 kg
Weight of case	4.47 kg
EAN of inner box	4038917086752
EAN case	4042836001574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-15

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**