




**Product Specification**  
**Snap-on lid, for 26 ml container 75.568, PE,  
transparent**

**Page 1**



**Product description**

Order number	76.569
Product description	Snap-on lid, for 26 ml container 75.568, transparent, material: PE, snap-on lid

**Product characteristics**

Closure type	snap-on lid
Suitable for	26 ml container 75.568

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Snap-on lid, for 26 ml container 75.568, PE,  
transparent**

**Page 2**

**Packaging**

Minimum order qty.	5000
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0016 kg
Weight of case	8.49 kg
EAN of inner box	4038917082556
EAN case	4038917001632

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**