



Product description

Order number	76.571
Product description	Snap-on lid, for 100-ml container 75.570, material: PS, transparent, collection and storage of urine, 100 piece(s)/tube bag

Product characteristics

Type of cap	Snap-on lid
Application	collection and storage of urine
Example of use	cap for urine collection tube
Suitable for	100-ml container 75.570

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	IVD product
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of inner packaging	tube bag
Piece(s) / inner box	1000
Piece(s) / case	1000
Piece(s) / pallet	32000
Case volume	0,0392 cbm
Weight of product	0,0009 kg
Weight of case	1,46 kg
EAN of inner box	4038917082570
EAN case	4038917001656

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

