



### Product description

Order number	51.9923.820
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, skirted conical base, transparent, cap assembled, with paper label, preparation: Boric acid, non-sterile, 500 piece(s)/case

### Product characteristics

Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	white/black
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, Boric acid, non-sterile**

**Page 2**

**Size**

Sample volume	25 ml
Diameter	25 mm
Length of product	90 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	red

**Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	non-sterile
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of inner packaging	case
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	15000
Case volume	0.059 cbm
Weight of product	0.009 kg
Weight of case	5.18 kg
EAN of inner box	4038917154550
EAN case	4038917030090

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-11-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**