



Product description

Order number	57.462
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	2500 x g

Size

Diameter	17 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Depth of case 398 mm
Width of case 398 mm
Height of case 310 mm
Case volume 0.0491 cbm
Weight of product 0.0063 kg
Weight of case 6.98 kg
EAN of inner box 4038917081702
EAN case 4038917000819

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

