



Product description

Order number	57.462.001
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, sterile, 100 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	2500 x g

Size

Diameter	17 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0077 kg
Weight of case	4.35 kg
EAN of inner box	4038917381048
EAN case	4038917381031

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

