



Product Specification
Urine Monovette® , 10 ml, cap blue, (LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 1



Product description

Order number	10.252.023
Product description	Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap blue, incl. collection tip, with plastic label, sterile, 1 piece(s)/blister

Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label
Colour of print/label	yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Example of use	Urine collection, transport, and storage
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



Product Specification
Urine Monovette® , 10 ml, cap blue, (LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 2

Recommended centrifugation 400 x g
Centrifugation time 5 Min
Connection Luer (m)

Size

Sample volume 10 ml
Diameter 15 mm
Length including cap 117 mm
Length excluding cap 102 mm
Length including cap and plunger 127 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap blue
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest blister
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 500
Piece(s) / pallet 13500
Depth of box 351 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Urine Monovette® , 10 ml, cap blue, (LxØ): 102 x 15 mm, 1 piece(s)/blister

Page 3

Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0084 kg
Weight of case	5.82 kg
EAN of inner box	4038917104456
EAN case	4038917009898

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT