



**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag**

**Page 1**



**Product description**

Order number	10.252.030
Product description	Urine Monovette® , nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, white, with light protection, graduated, material: PP, Luer screw cap, cap yellow, with plastic label, 64 piece(s)/bag

**Product characteristics**

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label
Colour of print/label	transparent/yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	yes
Application	Urine tube
Example of use	Urine collection, transport, and storage
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag**

**Page 2**

Recommended centrifugation 400 x g  
Centrifugation time 5 Min  
Connection Luer (m)

**Size**

Sample volume 10 ml  
Diameter 15 mm  
Length including cap 117 mm  
Length excluding cap 102 mm  
Length including cap and plunger 127 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product white  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Batched yes

**Packaging**

Minimum order qty. 512  
Type of smallest subpackaging bag  
Piece(s) / inner box 64  
Piece (s) / outer case 512  
Piece(s) / pallet 24576  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 205 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 10 ml, cap yellow, (LxØ): 102 x 15 mm, 64 piece(s)/bag**

Case volume	0.0301 cbm
Weight of product	0.0085 kg
Weight of case	4.85 kg
EAN of inner box	4038917279987
EAN case	4038917279840

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

