



Product Specification
Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag

Page 1



Product description

Order number	10.253.001
Product description	Urine Monovette® , Boric acid, nominal volume: 10 ml, Boric acid, Urine tube, (LxØ): 102 x 15 mm, transparent, material: PP, Luer screw cap, cap green, incl. collection tip, with paper label, 64 piece(s)/bag

Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Preparation concentration	15 mg/ml
Label/ Print	with paper label
Colour of print/label	green
Graduation	no
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base
Volume of work	10 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag**

Page 2

Connection Luer (m)

Size

Sample volume	10 ml
Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Safety data sheet	yes
Batched	yes

Packaging

Minimum order qty.	512
Type of smallest subpackaging	bag
Piece(s) / inner box	64
Piece (s) / outer case	512
Piece(s) / pallet	24576
Depth of case	388 mm
Width of case	378 mm
Height of case	205 mm
Case volume	0.0301 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Urine Monovette® , Boric acid, 10 ml, cap green,
(LxØ): 102 x 15 mm, 64 piece(s)/bag**

Page 3

Weight of product	0.0086 kg
Weight of case	4.9 kg
EAN of inner box	4038917280013
EAN case	4038917279871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT