



## Product Specification

**Urine Monovette® , Boric acid, 3.2 ml, cap green,  
(LxØ): 75 x 13 mm, 64 piece(s)/bag**

**Page 1**



### Product description

Order number	10.256.001
Product description	Urine Monovette® , Boric acid, nominal volume: 3.2 ml, Boric acid, Urine tube, (LxØ): 75 x 13 mm, transparent, fill mark, material: PP, Luer screw cap, cap green, incl. collection tip, with paper label, suitable for Tempus® sample transport system, 64 piece(s)/bag

### Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus® sample transport system
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Urine Monovette® , Boric acid, 3.2 ml, cap green,  
(LxØ): 75 x 13 mm, 64 piece(s)/bag**

**Page 2**

Volume of work	3,2 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

### Size

Sample volume	3.2 ml
Diameter	13 mm
Length including cap	90 mm
Length excluding cap	75 mm
Length including cap and plunger	109 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Safety data sheet	yes
Batched	yes

### Packaging

Minimum order qty.	512
Type of smallest subpackaging	bag
Piece(s) / inner box	64
Piece (s) / outer case	512
Piece(s) / pallet	36864
Depth of case	388 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , Boric acid, 3.2 ml, cap green,**  
**(LxØ): 75 x 13 mm, 64 piece(s)/bag**

**Page 3**

Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.1 kg
EAN of inner box	4038917982825
EAN case	4038917982832

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**