



## Product description

Order number	10.258.388
Product description	Urine Monovette® , nominal volume: 8.5 ml, without preparation, Urine tube, (LxØ): 92 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix, sterile, 1 piece(s)/blister

## Product characteristics

Type of collection	Aspiration
Type of preparation	without preparation
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/yellow
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 8.5 ml, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 2**

Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

**Size**

Sample volume	8.5 ml
Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	20000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine Monovette® , 8.5 ml, cap yellow, (LxØ): 92 x 15 mm, 1 piece(s)/blister**

Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0078 kg
Weight of case	5.5 kg
EAN of inner box	4038917949941
EAN case	4038917949958

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

