



## Product Specification

Urine Monovette® , Boric acid, 8.5 ml, cap green,  
(LxØ): 92 x 15 mm, 1 piece(s)/blister

Page 1



### Product description

Order number	10.260.021
Product description	Urine Monovette® , Boric acid, nominal volume: 8.5 ml, Boric acid, Urine tube, (LxØ): 92 x 15 mm, transparent, fill mark, Luer screw cap, cap green, incl. collection tip, with paper label, suitable for Tempus® sample transport system, sterile, 1 piece(s)/blister

### Product characteristics

Type of collection	Aspiration
Type of preparation	Boric acid
Label/ Print	with paper label
Colour of print/label	green/white
Graduation	fill mark
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	Urine tube
Suitable for	Tempus® sample transport system
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-12

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine Monovette® , Boric acid, 8.5 ml, cap green,  
(LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 2**

Volume of work	8,5 ml
Recommended centrifugation	400 x g
Centrifugation time	5 Min
Connection	Luer (m)

### Size

Sample volume	8.5 ml
Diameter	15 mm
Length including cap	106 mm
Length excluding cap	92 mm
Length including cap and plunger	127 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Safety data sheet	yes
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-12

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Urine Monovette® , Boric acid, 8.5 ml, cap green,  
(LxØ): 92 x 15 mm, 1 piece(s)/blister**

**Page 3**

Depth of box	350 mm
Width of box	150 mm
Height of box	170 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.008 kg
Weight of case	5.58 kg
EAN of inner box	4038917993685
EAN case	4038917993692

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-12

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**