



Product description

| | |
|---------------------|------------------------------------------------------------------------------------------|
| Order number | 11.1240.100 |
| Product description | Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag |

Product characteristics

| | |
|--------------------|-------------------------------|
| Type of collection | Vacuum |
| Light protection | no |
| Design | With wings on the tube holder |
| Application | Urine transfer |
| Suitable for | V-Monovette® Urine |
| Connection | Vacuum tube |

Size

| | |
|----------------|-------|
| Length of tube | 10 cm |
|----------------|-------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product yellow

Purity & certification

Product category In-vitro diagnostic (IVD)

CE certified CE supplier

Safety data sheet no

Batched yes

Packaging

Minimum order qty. 500

Type of smallest
subpackaging minigrip bag

Piece(s) / inner box 50

Piece (s) / outer case 500

Piece(s) / pallet 30000

Depth of box 300 mm

Width of box 250 mm

Height of box 10 mm

Depth of case 380 mm

Width of case 350 mm

Height of case 170 mm

Case volume 0.0226 cbm

Weight of product 0.002 kg

Weight of case 1.8 kg

EAN of inner box 6947462401157

EAN case 6947462401423

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

