



Product description

Order number	11.1240.100
Product description	Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag

Product characteristics

Type of collection	Vacuum
Light protection	no
Design	With wings on the tube holder
Application	Urine transfer
Suitable for	V-Monovette® Urine
Connection	Vacuum tube

Size

Length of tube	10 cm
----------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product yellow

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE supplier
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 500
Type of smallest minigrip bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 30000
Depth of box 300 mm
Width of box 250 mm
Height of box 10 mm
Depth of case 380 mm
Width of case 350 mm
Height of case 170 mm
Case volume 0.0226 cbm
Weight of product 0.002 kg
Weight of case 1.8 kg
EAN of inner box 6947462401157
EAN case 6947462401423

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

