



Product Specification
V-Monovette® Urine, 6 ml, cap yellow, (LxØ): 100
x 13 mm, 50 piece(s)/bag

Page 1



Product description

Order number	11.2352.001
Product description	V-Monovette® Urine, nominal volume: 6 ml, without preparation, vacuum system for closed urine transfer, (LxØ): 100 x 13 mm, transparent, fill mark, material: PET, membrane cap, cap yellow, with paper label, sterile, 50 piece(s)/bag

Product characteristics

Type of collection	Vacuum
Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	membrane cap
Closure type	push cap
Light protection	no
Application	vacuum system for closed urine transfer
Base shape	round base
Volume of work	6 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**V-Monovette® Urine, 6 ml, cap yellow, (LxØ): 100
x 13 mm, 50 piece(s)/bag**

Page 2

Recommended centrifugation 400 x g
Centrifugation time 5 Min
Connection V-Monovette®

Size

Sample volume 6 ml
Diameter 13 mm
Length including cap 106 mm
Length excluding cap 100 mm

Material & colours

Product material Polyethylene therephthalate (PET)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow
Membrane material rubber
Tube material Polyethylene therephthalate (PET)

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Purity standard sterile
Sterilisation Gamma irradiation
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 34000
Depth of case 428 mm
Width of case 398 mm
Height of case 135 mm
Case volume 0.023 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
V-Monovette® Urine, 6 ml, cap yellow, (LxØ): 100
x 13 mm, 50 piece(s)/bag

Page 3

Weight of product	0.0072 kg
Weight of case	4.14 kg
EAN of inner box	4038917245098
EAN case	4038917245081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT