



### **Product description**

Order number	14.1217
Product description	Haemo-Diff® , without smear edge, for producing blood smears, 500 piece(s)/case

### **Product characteristics**

Design	without smear edge
Application	Blood smear

### **Size**

Length of product	26,5 mm
-------------------	---------

### **Material & colours**

Product material	Styrene-butadiene (SB)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of inner packaging	case
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	648000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0,0051 cbm
Weight of product	0,0003 kg
Weight of case	0,81 kg
EAN of inner box	4038917085823
EAN case	4038917004459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

