



### **Product description**

Order number	14.1217
Product description	Haemo-Diff® , without smear edge, for producing blood smears, 500 piece(s)/case

### **Product characteristics**

Design	without smear edge
Application	Blood smear

### **Size**

Length of product	26,5 mm
-------------------	---------

### **Material & colours**

Product material	Styrene-butadiene (SB)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category                    no medical device | IVD  
Batched                                yes

### **Packaging**

Minimum order qty.                2000  
Type of inner packaging            case  
Piece(s) / inner box                500  
Piece(s) / case                        2000  
Piece(s) / pallet                    648000  
Depth of case                        268 mm  
Width of case                        191 mm  
Height of case                        100 mm  
Case volume                         0,0051 cbm  
Weight of product                  0,0003 kg  
Weight of case                        0,81 kg  
EAN of inner box                    4038917085823  
EAN case                                4038917004459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

