



Product description

Order number	14.1218
Product description	Haemo-Diff® , with smear edge, for producing blood smears, 500 piece(s)/case

Product characteristics

Design	with smear edge
Application	Blood smear
Connection	S-Monovette®

Size

Length of product	21 mm
-------------------	-------

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2000
Type of inner packaging case
Piece(s) / inner box 500
Piece(s) / case 2000
Piece(s) / pallet 216000
Depth of box 191 mm
Width of box 129 mm
Height of box 129 mm
Depth of case 398 mm
Width of case 268 mm
Height of case 150 mm
Case volume 0,016 cbm
Weight of product 0,0013 kg
Weight of case 3,27 kg
EAN of inner box 4038917085830
EAN case 4038917004466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

