




**Product Specification**  
**Haemo-Diff® , with smear edge, for producing**  
**blood smears**

**Page 1**



**Product description**

Order number	14.1218
Product description	Haemo-Diff® , with smear edge, for producing blood smears, 500 piece(s)/case

**Product characteristics**

Design	with smear edge
Application	Blood smear
Connection	S-Monovette®

**Size**

Length of product	21 mm
-------------------	-------

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Haemo-Diff® , with smear edge, for producing**  
**blood smears**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	216000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	398 mm
Width of case	268 mm
Height of case	150 mm
Case volume	0.016 cbm
Weight of product	0.0013 kg
Weight of case	3.27 kg
EAN of inner box	4038917085830
EAN case	4038917004466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-23

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**