



Product description

Order number	14.1513.002
Product description	Colour-coded cap, for S-Monovette® , yellow, 100 piece(s)/bag

Product characteristics

Suitable for	S-Monovette®
--------------	--------------

Size

Length of product	13 mm
-------------------	-------

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	yellow

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	324000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0003 kg
Weight of case	0.47 kg
EAN of inner box	4038917251556
EAN case	4038917251549

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

