




Product Specification
Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation

Page 1



Product description

Order number	53.419
Product description	Seraplas® valve filter, V13, blue, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

Product characteristics

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

Size

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Seraplas® valve filter, blue, for separation of
serum/plasma from the blood cells after
centrifugation**

Page 2

Material & colours

Product material	Polyethylene (PE)
Colour of product	blue

Purity & certification

Product category	no medical device IVD
Safety data sheet	no
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	75000
Depth of box	119 mm
Width of box	99 mm
Height of box	98 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0016 kg
Weight of case	2.15 kg
EAN of inner box	4038917081481
EAN case	4038917000574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT