



**Product Specification**  
**Seraplas® valve filter, yellow, for separation of  
serum/plasma from the blood cells after  
centrifugation**

**Page 1**



**Product description**

Order number	53.428
Product description	Seraplas® valve filter, V15, yellow, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

**Product characteristics**

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

**Size**

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Seraplas® valve filter, yellow, for separation of serum/plasma from the blood cells after centrifugation**

**Page 2**

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	yellow

**Purity & certification**

Product category	no medical device   IVD
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	75000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.002 kg
Weight of case	2.61 kg
EAN of inner box	4038917081528
EAN case	4038917000611

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**