




**Product Specification**  
**Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation**

**Page 1**



**Product description**

Order number	53.677
Product description	Seraplas® valve filter, V11, white, for separation of serum/plasma from the blood cells after centrifugation, 500 piece(s)/bag

**Product characteristics**

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

**Size**

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation**

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	white

**Purity & certification**

Product category	no medical device   IVD
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0015 kg
Weight of case	8.5 kg
EAN of inner box	4038917083317
EAN case	4038917002257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

