




**Product Specification**  
**Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation**



**Product description**

Order number 53.677  
Product description Seraplas® valve filter, V11, white, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/bag

**Product characteristics**

Application venous blood collection  
Suitable for separation of serum/plasma from the blood cells after centrifugation

**Size**

Length of product 92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Seraplas® valve filter, white, for separation of serum/plasma from the blood cells after centrifugation**

**Material & colours**

Product material Polyethylene (PE)  
Colour of product white

**Purity & certification**

Product category no medical device | IVD  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 5000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 5000  
Piece(s) / pallet 120000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 415 mm  
Case volume 0.0657 cbm  
Weight of product 0.0015 kg  
Weight of case 8.5 kg  
EAN of inner box 4038917083317  
EAN case 4038917002257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

