



Product description

Order number	53.9924.613
Product description	Serum filter, translucent, (LxØ): 100 x 16 mm, filtering fibrin from serum/plasma of centrifuged blood, 100 piece(s)/bag

Product characteristics

Application	Serum filter
Example of use	filtering fibrin from serum/plasma of centrifuged blood

Size

Diameter	16 mm
Length of product	100 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	translucent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	60000
Depth of case	428 mm
Width of case	358 mm
Height of case	234 mm
Case volume	0.0357 cbm
Weight of case	3 kg
EAN of inner box	4038917155441
EAN case	4038917398336

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

