



Product description

Order number	55.1570.002
Product description	Adapter tube, (LxØ): 54 x 11 mm, material: PP, brown, 250 piece(s)/bag

Product characteristics

Light protection	yes
Suitable for	S-Monovette® with Ø 8 mm
Base shape	flat base

Size

Diameter	11 mm
Length of product	54 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Adapter tube, (LxØ): 54 x 11 mm, PP, brown

Page 2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0012 kg
Weight of case	1.54 kg
EAN of inner box	4038917131018
EAN case	4038917023573

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT