



### Product description

Order number	55.1570
Product description	Adapter tube, (LxØ): 54 x 11 mm, material: PP, transparent, 250 piece(s)/bag

### Product characteristics

Suitable for	S-Monovette® with Ø 8 mm
Base shape	flat base

### Size

Diameter	11 mm
Length of product	54 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	78000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0.018 cbm
Weight of product	0.0011 kg
Weight of case	1.51 kg
EAN of inner box	4038917087063
EAN case	4038917005531

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

