



### **Product description**

Order number	55.1571
Product description	Adapter tube, (LxØ): 65 x 13 mm, material: PP, transparent, 250 piece(s)/bag

### **Product characteristics**

Closure type	screw cap
Suitable for	S-Monovette® with Ø 8 mm, Microvette® and micro tubes
Base shape	round base

### **Size**

Diameter	13 mm
Length of product	65 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                    no

**Packaging**

Minimum order qty.                      1000  
Type of smallest                            bag  
subpackaging  
Piece(s) / inner box                      250  
Piece (s) / outer case                    1000  
Piece(s) / pallet                          72000  
Depth of case                              333 mm  
Width of case                               303 mm  
Height of case                              190 mm  
Case volume                                0.0192 cbm  
Weight of product                        0.0019 kg  
Weight of case                              2.23 kg  
EAN of inner box                          4038917242523  
EAN case                                    4038917242516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

