



Product description

Order number	55.1572
Product description	Adapter tube, (LxØ): 55 x 13 mm, material: PP, transparent, 250 piece(s)/bag

Product characteristics

Closure type	screw cap
Suitable for	S-Monovette® with Ø 8 mm
Base shape	round base

Size

Diameter	13 mm
Length of product	55 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	72000
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0014 kg
Weight of case	1.71 kg
EAN of inner box	4038917283823
EAN case	4038917283816

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

