



### Product description

Order number	55.1577
Product description	Carrier tube, (LxØ): 60 x 11.5 mm, material: PP, transparent, 250 piece(s)/bag

### Product characteristics

Type of cap	screw cap
Suitable for	carrier tube for Multivette®

### Size

Diameter	11,5 mm
Length of product	60 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Purity & certification**

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	78000
Depth of case	428 mm
Width of case	323 mm
Height of case	130 mm
Case volume	0,018 cbm
Weight of product	0,0013 kg
Weight of case	1,7 kg
EAN of inner box	4038917087070
EAN case	4038917005548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

