



Product description

Order number	55.525.003
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, suitable for Microvette® and micro tube with screw cap, 250 piece(s)/bag

Product characteristics

Suitable for	Microvette® and micro tube with screw cap
Base shape	round base
Volume of work	5 ml

Size

Diameter	13 mm
Length of product	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PP

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest
subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 2000
Piece(s) / pallet 72000
Depth of case 398 mm
Width of case 398 mm
Height of case 265 mm
Case volume 0.042 cbm
Weight of product 0.0026 kg
Weight of case 5.7 kg
EAN of inner box 4038917109154
EAN case 4038917012362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT