



### **Product description**

Order number	60.503
Product description	Screw cap tube, working volume: 6.5 ml, (LxØ): 90 x 13 mm, material: PP, round base, transparent, screw cap, natural, cap enclosed, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Light protection	no
Base shape	round base
Volume of work	6,5 ml
Centrifugation max (RCF)	3000 x g

### **Size**

Diameter	13 mm
Height of product	90 mm
Length of product	90 mm
Length including cap	106 mm
Length excluding cap	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Safety data sheet	no
Batched	no

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	458 mm
Width of case	408 mm
Height of case	260 mm
Case volume	0.0486 cbm
Weight of product	0.0034 kg
Weight of case	7.46 kg
EAN of inner box	4038917555838
EAN case	4038917555821

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

