



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
without cap, 100 piece(s)/bag**

**Page 1**



### Product description

Order number	60.612.010
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, without cap, 100 piece(s)/bag, 1,000 piece(s)/carton

### Product characteristics

Graduation	no
Base shape	conical false bottom, rounded tube bottom
Volume of work	5 ml

### Size

Diameter	15.3 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
without cap, 100 piece(s)/bag**

**Page 2**

### Material & colours

Product material Polypropylene (PP)  
Colour of product transparent

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### Packaging

Minimum order qty. 1000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 36000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 265 mm  
Case volume 0.042 cbm  
Weight of product 0.004 kg  
Weight of case 4.54 kg  
EAN of inner box 4038917114936  
EAN case 4038917014458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**