



Product Specification

**Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 1,000 piece(s)/bag**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 60.614.010 |
| Product description | Screw cap tube, working volume: 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton |

Product characteristics

| | |
|----------------|---|
| Graduation | yes |
| Base shape | conical false bottom, rounded tube bottom |
| Volume of work | 2,5 ml |

Size

| | |
|-------------------|-------|
| Diameter | 13 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-01-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Screw cap tube, 2.5 ml, (LxØ): 75 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 1,000 piece(s)/bag**

Page 2

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | In-vitro diagnostic, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 458 mm |
| Width of case | 408 mm |
| Height of case | 260 mm |
| Case volume | 0.0486 cbm |
| Weight of product | 0.0024 kg |
| Weight of case | 5.45 kg |
| EAN of inner box | 4038917353182 |
| EAN case | 4038917353175 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-01-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT