



Product Specification

**Screw cap tube, 2.5 ml, 2.5 ml, (LxØ): 75 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 500 piece(s)/StackPack**

Page 1



Product description

Order number	60.614.015
Product description	Screw cap tube, working volume: 2.5 ml, nominal volume: 2.5 ml, (LxØ): 75 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, without cap, 500 piece(s)/StackPack, 2,000 piece(s)/carton

Product characteristics

Graduation	yes
Base shape	conical false bottom, rounded tube bottom
Volume of work	2,5 ml

Size

Sample volume	2.5 ml
Diameter	13 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-02-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Screw cap tube, 2.5 ml, 2.5 ml, (LxØ): 75 x 13 mm,
conical false bottom, rounded tube bottom, PP,
without cap, 500 piece(s)/StackPack**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	112000
Depth of box	413 mm
Width of box	76 mm
Height of box	189 mm
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.0024 kg
Weight of case	5.57 kg
EAN of inner box	4038917280044
EAN case	4038917280037

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-02-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT