



Product description

Order number	62.503.210
Product description	Screw cap tube, working volume: 6.5 ml, Ø: 13 mm, round base, transparent, without cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Light protection	no
Base shape	round base
Volume of work	6,5 ml
Centrifugation max (RCF)	3000 x g

Size

Diameter	13 mm
Height of product	90 mm

Material & colours

Colour of product	transparent
Tube material	Polypropylene (PP)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	50000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap tube, 6.5 ml, Ø: 13 mm, with print

Page 2

Case volume	0.0475 cbm
Weight of product	0.0027 kg
Weight of case	6.12 kg
EAN of inner box	4038917429061
EAN case	4038917429054

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-12-10

This document was prepared by EDP support and is valid without signature.



SARSTEDT