



Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, rounded tube bottom, PP,
cap assembled, 100 piece(s)/bag**

Page 1



Product description

Order number	62.617
Product description	Screw cap tube, working volume: 3.5 ml, (LxØ): 92 x 13 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with graduations, cap assembled, 100 piece(s)/bag, 1,000 piece(s)/carton

Product characteristics

Graduation	yes
Base shape	conical false bottom, rounded tube bottom
Volume of work	3,5 ml

Size

Diameter	13 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Screw cap tube, 3.5 ml, (LxØ): 92 x 13 mm,
conical false bottom, rounded tube bottom, PP,
cap assembled, 100 piece(s)/bag**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Depth of case 458 mm
Width of case 408 mm
Height of case 205 mm
Case volume 0.0383 cbm
Weight of product 0.0036 kg
Weight of case 4.34 kg
EAN of inner box 4038917082921
EAN case 4038917079150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT