



Product description

Order number 65.1379.002
Product description Cap, red, suitable for tubes Ø 12-17 mm, 1,000 piece(s)/bag

Product characteristics

Cap standard cap
Suitable for tubes Ø 12-17 mm

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product red

Purity & certification

Product category no medical device | IVD
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	240000
Case volume	0.0441 cbm
Weight of product	0.0004 kg
Weight of case	5.36 kg
EAN of inner box	4038917127981
EAN case	4038917021517

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

