




Product Specification
Archiving cap, light blue, suitable for
S-Monovette® , tubes Ø 13-16 mm

Page 1



Product description

Order number	65.647.020
Product description	Archiving cap, with filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag

Product characteristics

Suitable for	S-Monovette® , tubes Ø 13-16 mm
Filter	yes

Material & colours

Product material	Polypropylene (PP)
Colour of product	light blue
Colour of cap	light blue

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Archiving cap, light blue, suitable for
S-Monovette® , tubes Ø 13-16 mm

Page 2

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.001 kg
Weight of case	5.37 kg
EAN of inner box	4038917250085
EAN case	4038917250078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT