



**Product Specification**  
**Archiving cap, light blue, suitable for**  
**S-Monovette® , tubes Ø 13-16 mm**

**Page 1**



**Product description**

Order number	65.647
Product description	Archiving cap, without filter, light blue, suitable for S-Monovette® , tubes Ø 13-16 mm, 1,000 piece(s)/bag

**Product characteristics**

Suitable for	S-Monovette® , tubes Ø 13-16 mm
Filter	no

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	light blue

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-23

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Archiving cap, light blue, suitable for**  
**S-Monovette<sup>®</sup>, tubes Ø 13-16 mm**

**Page 2**

**Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0009 kg
Weight of case	5.2 kg
EAN of inner box	4038917158640
EAN case	4038917313292

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-23

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**