



Product description

| | |
|---------------------|--|
| Order number | 65.803 |
| Product description | Push cap, natural, suitable for tubes Ø 15.7 mm & cuvette 67.743, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------|----------------------------------|
| Cap | standard cap |
| Closure type | push cap |
| Suitable for | tubes Ø 15.7 mm & cuvette 67.743 |

Size

| | |
|-------------------|----------|
| Diameter | 16.35 mm |
| Height of product | 20.45 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Low Density Polyethylene (LD-PE) |
| Colour of product | natural |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 5000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 5000
Piece(s) / pallet 210000
Depth of case 398 mm
Width of case 398 mm
Height of case 250 mm
Case volume 0.0396 cbm
Weight of product 0.0013 kg
Weight of case 7.1 kg
EAN of inner box 4038917084192
EAN case 4038917003056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

