



Product Specification
Insert, suitable for tubes and S-Monovette® Ø 16 mm, transparent

Page 1



Product description

Order number 73.1053
Product description Insert, suitable for tubes and S-Monovette® Ø 16 mm, material: PS, length: 23.5 mm, transparent, 500 piece(s)/bag

Product characteristics

Graduation no
Suitable for tubes and S-Monovette® Ø 16 mm
Base shape skirted conical base

Size

Diameter of opening 11 mm
Length of product 23.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Insert, suitable for tubes and S-Monovette® Ø 16 mm, transparent

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 156000
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm
Case volume 0.018 cbm
Weight of product 0.0009 kg
Weight of case 2.26 kg
EAN of inner box 4038917085168
EAN case 4038917003841

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-01-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT