




Product Specification
Wall holder, for tournistrip® disposable
tourniquet (item no 95.1006)

Page 1



Product description

Order number 95.1005
Product description Wall holder, for tournistrip® disposable tourniquet (item no 95.1006)

Product characteristics

Suitable for tournistrip® disposable tourniquet (item no 95.1006)

Material & colours

Product material metal

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Wall holder, for tournistrip® disposable
tourniquet (item no 95.1006)

Page 2

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.21 kg
Weight of case	0.33 kg
EAN of inner box	4038917348560
EAN case	4038917348553

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT