



Product description

Order number	95.1006
Product description	tournistrip® disposable tourniquet, material: special paper/plastic composite material, latex-free, 200 piece(s)/dispenser

Product characteristics

Application	venous blood collection
-------------	-------------------------

Material & colours

Product material	special paper/plastic composite material
------------------	--

Purity & certification

Product category	Medical device
CE certified	CE supplier
Latex-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	200
Type of inner packaging	dispenser
Piece(s) / inner box	200
Piece(s) / case	4800
Piece(s) / pallet	201600
Depth of box	181 mm
Width of box	36 mm
Height of box	181 mm
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0023 kg
Weight of case	12.43 kg
EAN of inner box	5060728810157
EAN case	5060728810164

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

