



Product description

Order number	95.1007.001
Product description	Tourniquet, material: PL, latex-free, 1 piece(s)/bag

Product characteristics

Colour of print/label	white
Application	venous blood collection

Material & colours

Product material	polyester (PL)
Colour of product	red

Purity & certification

Product category	Medical device
CE certified	CE supplier
Latex-free	yes
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	200
Piece(s) / pallet	12800
Depth of case	330 mm
Width of case	330 mm
Height of case	220 mm
Case volume	0.024 cbm
Weight of product	0.0306 kg
Weight of case	6.86 kg
EAN of inner box	4260150884707
EAN case	4260150884806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

