



### **Product description**

Order number 95.1007.001  
Product description Tourniquet, material: PE, latex-free, 1 piece(s)/bag

### **Product characteristics**

Colour of print/label white  
Application venous blood collection

### **Material & colours**

Product material polyester (PE)  
Colour of product red

### **Purity & certification**

Product category Medical device  
CE certified CE supplier  
Latex-free yes  
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Packaging**

Minimum order qty.	1
Type of inner packaging	bag
Piece(s) / inner box	1
Piece(s) / case	200
Piece(s) / pallet	12800
Case volume	0.0246 cbm
Weight of product	0.0306 kg
Weight of case	6.86 kg
EAN of inner box	4260150884707
EAN case	4260150884806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

