



Product description

Order number 95.1007.101
Product description Tourniquet, latex-free, 250 piece(s)/bag

Product characteristics

Application venous blood collection

Material & colours

Colour of product yellow

Purity & certification

Product category no medical device | IVD
Latex-free yes
Batched yes

Packaging

Minimum order qty. 250
Type of smallest subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 2000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tourniquet, latex-free

Page 2

Piece(s) / pallet	42000
Case volume	0 cbm
Weight of case	7.01 kg
EAN of inner box	4038917263214
EAN case	4038917398220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT