



Product description

Order number 95.1007.101
Product description Tourniquet, latex-free

Product characteristics

Application venous blood collection

Purity & certification

Product category no medical device | IVD
Latex-free yes
Batched yes

Packaging

Minimum order qty. 250
Type of smallest case
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 2000
Piece(s) / pallet 42000
Case volume 0 cbm
Weight of case 7.01 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification
Tourniquet, latex-free**

Page 2

EAN of inner box 4038917263214
EAN case 4038917398220

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT